


HUNTINGTON DISEASE KIT 1 - FL

System for genetic testing for Huntington Disease (HD)
by tripletrepeat primed PCR (TP-PCR) and
fragment sizing through capillary electrophoresis



CytoGen GmbH
Langgasse 73
35576 Wetzlar, Germany

Tel. + 49 6441-6795588
Fax +49 6441-6795589

web: www.cytogen.info
e-mail: cytogen@eurobiz.de

Huntington Disease is a progressive neurodegenerative disorder that presents with motor symptoms, cognitive impairment and psychiatric disturbances. The first symptoms usually manifest between 35 and 50 years of age and duration of the disease is between 15 and 20 years. A small number of cases present before the age of 20 (juvenile onset) and about 25% of cases present after 50 years of age.

It is characterized pathologically by loss of specific neuronal populations in many brain regions, although the pathology is not limited to neurons. Neuropathological features include selective degeneration of neurons in the caudate and putamen and less severe loss in the cerebral cortex.

HD is caused by the expansion of an unstable polymorphic **trinucleotide (CAG)_n** repeat in exon 1 of the **HTT gene** (4p16.3), which translates into an extended polyglutamine tract in the protein.

Alleles with <27 CAG repeats are classified as **normal**, whereas alleles with **≥36 repeats** are detected in patients. Alleles with **27-35** repeats (called **mutable** or intermediated alleles) are not associated with the disease symptoms but can expand into the affected range upon (predominantly paternal) germline transmission and thus cause HD in offspring.

Repeats of **36-39 CAG** are **incompletely penetrant** and can be found in affected individual as well as individuals who show no clinical symptoms. HD appears to be fully penetrant for allele size of **≥ 40** repeats.

Genetic testing for HD is usually requested by a neurologist or a clinical geneticist for a confirmation or exclusion of a clinical diagnosis; if the diagnosis is confirmed by DNA analysis, the patient and family members should be referred for genetic counseling and a possible offer of a presymptomatic testing.

Containing of the kit

Label	Contents
D4/6FAM-HD	Mix for the amplification of the FMR-1 gene
MASTER MIX	Mix for the amplification of the control gene
Mix for amplification of trinucleotide	Wild-type Control

How does the kit work?

The **Huntington Disease kit-FL** is used for the molecular testing of the Huntington disease. This system amplifies the **HD (CAG)_n** repeat region by TP-PCR (triplet-repeat primed PCR) using a fluorescently labeled forward primer located upstream of the (CAG)_n region and a chimeric reverse primer located partially within the (CAG)_n region.

The chimeric reverse primer hybridizes to multiple locations within the (CAG)_n repeat region, creating PCR products of varying sizes. Reactions are separated by capillary electrophoresis. TP-PCR provides a characteristic ladder on the fluorescence trace, enabling the rapid identification of a large pathogenic repeats that cannot be amplified using flanking primers. The **»HD« sequence » specific » standard** (provided), that has to be included on each run, allows to detect accurately normal alleles, mutable alleles and HD alleles.

Starting samples: peripheral blood. **DNA isolation method:** QIAamp DNA blood mini kit, QIAcube, QIA-Asymphony (Qiagen), High Pure PCR template preparation kit (Roche). **DNA Sequencer:** CEQ 8000/8800, GeXP Genetic Analysis System (Beckman Coulter); 310, 3100, 3130, 3730, 3500 Genetic Analyzers (Applied Biosystems).

Procedure: according to "EMQN/CMGS best practice guidelines for molecular genetic testing of HD" 2013.

Product	Unit	Cat.-No.
Huntington disease Kit 1-FL	40 tests	HD.01FL

Distributor:

CytoGen GmbH
Langgasse 73
35576 Wetzlar, Germany
Tel. + 49 6441 6795588
Fax +49 6441 6795589
web: www.cytogen.info
e-mail: cytogen@eurobiz.de

Producer:

EXPERTEAM S.r.l.
via della Libertà 12
30175 Marghera, Italy
tel.: +39 041 5093101
fax: +39 041 5093102
web: www.experteam.it
e-mail: expertm@vegapark.ve.it